

DuoFLEX® MRI Coil Suite

VALUE FOR TECHNOLOGISTS

The DuoFLEX® MRI Coil Suite contours to your patient and your budget, offering flexibility without compromise. The suite includes two 10cm and two 24cm coils which provide versatility for the technician and comfort for the patient. The DuoFLEX design helps deliver high-resolution imaging, even on old scanners, that radiologists can depend on to make accurate diagnoses.

DRIVING EFFICIENCY & QUALITY PATIENT CARE

Getting it done right the first time creates a win-win experience for everyone. Managing day-to-day scheduling and throughput while ensuring patient comfort are your highest priorities. Spend less time between scans replacing and preparing different surface coils and produce more high-quality, diagnostic images. The DuoFLEX® MRI Coil Suite is flexible enough to scan any body part; if you can imagine it, you can scan it.

Benefits:

- **Patient comfort:** extra-long cables allow you to bring the coil to the patient versus moving the patient to the coil
- **Flex-fit paddles** contour to small and large anatomy which increases SNR
- **High resolution** for small-parts imaging with 10cm paddles, supports more accurate diagnoses of hard-to-scan anatomy
- **DuoFLEX® 24cm and 10cm flex paddles** work interchangeably for the best image and comfort outcomes

"These are the best coils I have ever used in terms of image quality, versatility and speed of scanning. DuoFLEX coils are capable of scanning most anything that will fit into their field of view! In addition to our musculoskeletal needs, we use them for female pelvis, hips, brains, brachial plexus, TMJ exams and much more!"

*- Grant Thompson
MRI Technologist*

